



## 3D flow focusing chips

### Product Datasheet

#### Contents

1	Description.....	4
2	Benefits.....	7
3	Specification.....	9
4	Geometry.....	11
4.1	3D Flow Focusing Chip - 65 $\mu\text{m}$ (P.N: 3200958, 3200959, 3200960).....	11
4.2	3D Flow Focusing Chip - 100 $\mu\text{m}$ (P.N: 3200433, 3200434, 3200515).....	11
4.3	3D Flow Focusing Chip - 170 $\mu\text{m}$ (P.N: 3200835, 3200836, 3200837).....	15
4.4	3D Flow Focusing Chip - 100 $\mu\text{m}$ - 2 Reagent (P.N: 3200435, 3200436).....	15
4.5	3D Flow Focusing Chip – 30 $\mu\text{m}$ (P.N.: 3200934, 3200935, 3200936).....	19
4.6	3D Flow Focusing Chip - 14 $\mu\text{m}$ (P.N: 3200437, 3200438).....	19
5	Surface Coatings .....	22
6	Droplet Formation .....	22
7	Custom Options .....	26
8	IP License .....	27



**Dolomite Microfluidics is a brand of Blacktrace**

Registered office: The Dolomite Centre Ltd. 1, Anglian Business Park, Royston, SG8 5TW, UK

Reg in England No. 04257809





Part Name	Part Number
3D Flow Focusing Chip – 65 $\mu\text{m}$ - Hydrophilic	3200958
3D Flow Focusing Chip – 65 $\mu\text{m}$ - Hydrophobic	3200959
3D Flow Focusing Chip – 65 $\mu\text{m}$ - Fluorophilic	3200960
3D Flow Focusing Chip – 100 $\mu\text{m}$ - Hydrophilic	3200433
3D Flow Focusing Chip – 100 $\mu\text{m}$ - Hydrophobic	3200434
3D Flow Focusing Chip – 100 $\mu\text{m}$ - Fluorophilic	3200515
3D Flow Focusing Chip – 170 $\mu\text{m}$ - Hydrophilic	3200835
3D Flow Focusing Chip – 170 $\mu\text{m}$ - Hydrophobic	3200836
3D Flow Focusing Chip – 170 $\mu\text{m}$ - Fluorophilic	3200837
3D Flow Focusing Chip – 100 $\mu\text{m}$ - 2 Reagent - Hydrophilic	3200435
3D Flow Focusing Chip – 100 $\mu\text{m}$ - 2 Reagent - Hydrophobic	3200436
3D Flow Focusing Chip – 30 $\mu\text{m}$ - Hydrophilic	3200934
3D Flow Focusing Chip – 30 $\mu\text{m}$ – Hydrophobic	3200935
3D Flow Focusing Chip – 30 $\mu\text{m}$ – Fluorophilic	3200936
3D Flow Focusing Chip – 14 $\mu\text{m}$ - Hydrophilic	3200437
3D Flow Focusing Chip – 14 $\mu\text{m}$ – Hydrophobic	3200438



**Dolomite Microfluidics is a brand of Blacktrace**

Registered office: The Dolomite Centre Ltd. 1, Anglian Business Park, Royston, SG8 5TW, UK

Reg in England No. 04257809





# 1 Description


3D flow focusing chips are glass microfluidic devices designed for generating droplets and bubbles. They are produced by wet etch lithography and feature channels that are 14  $\mu\text{m}$ , 30  $\mu\text{m}$ , 100  $\mu\text{m}$  and 170  $\mu\text{m}$  in etch depth. The precision of fabrication allows for creation of extremely monodisperse droplets and bubbles.

Certain fluids such as polymer solutions can cause surface fouling in standard flow focusing chips, reducing production run duration and chip lifetime. 3D flow focusing chips reduce this fouling by utilising a “pore” structure on the outlet side of the droplet-forming junction. This induces the droplet fluid to detach from the top and bottom walls well upstream of the breakoff point. As the length of the interface in contact with these walls is reduced, the rate of surface fouling is reduced, allowing much longer production runs.



*Figure 1 X-junction of 3D Flow Focusing Chip - 100  $\mu\text{m}$*


This range includes the following chips:



**3D Flow Focusing Chip - 65 µm:** suitable for generating droplets around 65 µm depending on fluid properties. Available in Hydrophilic (Part No. 3200958), Hydrophobic (Part No. 3200959) and Fluorophilic (Part no. 3200960) versions.

**3D Flow Focusing Chip - 100 µm and 170 µm:** suitable for generating droplets between 25-70 µm and 50-140 µm depending on fluid properties. Available in Hydrophilic (Part No. 3200433 and 3200835), Hydrophobic (Part No. 3200434 and 3200836) and Fluorophilic (Part No. 3200515 and 3200837) versions for generating organic-in-aqueous, aqueous-in-organic oil and aqueous-in-fluorocarbon oil emulsions respectively.

**3D Flow Focusing Chip - 100 µm - 2 Reagent:** similar to the 100 µm chip above, this chip features two separate inputs for the droplet fluid, allowing in-droplet mixing, in-droplet reactions or formation of Janus particles. Available in Hydrophilic (Part No. 3200435) and Hydrophobic (Part No. 3200436) versions.



**3D Flow Focusing Chip - 30 µm:** suitable for generating droplets between 10-28 µm depending on fluid properties. Available in Hydrophilic (Part No. 3200934), Hydrophobic (Part No. 3200935) and Fluorophilic (Part no. 3200936) versions.



**3D Flow Focusing Chip - 14 µm:** suitable for generating droplets between 5-12 µm depending on fluid properties. Available in Hydrophilic (Part No. 3200437) and Hydrophobic (Part No. 3200438) versions.

***N.B.** It may be possible to generate droplets of larger diameter than stated above. However, droplets with diameter near to or greater than the channel depth are likely to make contact the channel walls, leading to surface fouling in the outlet channel.*

*The 100 µm and 170 µm depth variants are compatible with our H Interface (Part No. 3000155) and Linear Connector 4-Way (Part No. 3000024). The 14µm and the 30µm depth variants are compatible with our Top Interface 4-Way (4mm) (Part No. 3000109) and Linear Connector 4-Way (Part No. 3000024).*

## 2 Benefits

- Extended lifetime when working with polymerising or surface fouling droplet fluids
- Quick connection using Dolomite Interfaces and Connectors
- Simple to use
- High visibility
- Good accessibility for microscope objective
- Extremely consistent droplet size
- Wide temperature and pressure range

- Excellent chemical compatibility
- Hydrophilic, Hydrophobic and Fluorophilic versions available



### 3 Specification

Specification	170 $\mu\text{m}$ 3D Flow Focusing Chip		100 $\mu\text{m}$ 3D Flow Focusing Chip		2 Reagent 100 $\mu\text{m}$ 3D Flow Focusing Chip		30 $\mu\text{m}$ 3D Flow Focusing Chip		65 $\mu\text{m}$ 3D Flow Focusing Chip		14 $\mu\text{m}$ 3D Flow Focusing Chip	
	3200835	3200836 3200837	3200433	3200434 3200515	3200435	3200436	3200934	3200935 3200936	3200958	3200959 3200960	3200437	3200438
Channel cross-section at junction (depth x width)	170 $\mu\text{m}$ x 175 $\mu\text{m}$		100 $\mu\text{m}$ x 105 $\mu\text{m}$				30 $\mu\text{m}$ x 33 $\mu\text{m}$		65 $\mu\text{m}$ x 70 $\mu\text{m}$		14 $\mu\text{m}$ x 17 $\mu\text{m}$	
Channel cross-section elsewhere (depth x width)	170 $\mu\text{m}$ x 175 $\mu\text{m}$		100 $\mu\text{m}$ x 105 $\mu\text{m}$				30 $\mu\text{m}$ x 1500 $\mu\text{m}$ (after junction)		65 $\mu\text{m}$ x 1500 $\mu\text{m}$		14 $\mu\text{m}$ x 500 $\mu\text{m}$	
Nominal pore dimensions (depth x width)	103 $\mu\text{m}$ x 140 $\mu\text{m}$		54 $\mu\text{m}$ x 75 $\mu\text{m}$				23 $\mu\text{m}$ x 16 $\mu\text{m}$		34 $\mu\text{m}$ x 49 $\mu\text{m}$		8 $\mu\text{m}$ x 12 $\mu\text{m}$	
Number of droplet phase inputs	1				2		1		1		1	
On-chip filter pore size (depth x width)	-				-		-		-		7 $\mu\text{m}$ x 17 $\mu\text{m}$	
Channel length after junction	15.4 mm				11.5 mm		24 mm		24.6 mm		27.4 mm	
Back pressure with 100 $\mu\text{L}/\text{min}$ flow (water)	100 mbar		121 mbar		116 mbar		-		150 mBar		-	
Back pressure with 10 $\mu\text{L}/\text{min}$ flow (water)	-				-		280 mbar		-		640 mbar	
Surface roughness of channels (Ra)					5nm							

Specification	170 $\mu$ m 3D Flow Focusing Chip		100 $\mu$ m 3D Flow Focusing Chip		2 Reagent 100 $\mu$ m 3D Flow Focusing Chip		30 $\mu$ m 3D Flow Focusing Chip		65 $\mu$ m 3D Flow Focusing Chip		14 $\mu$ m 3D Flow Focusing Chip	
	3200835	3200836 3200837	3200433	3200434 3200515	3200435	3200436	3200934	3200935 3200934	3200958	3200959 3200960	3200437	3200438
Chip size (length x width x thickness)	22.5mm x 15.0mm x 4mm											
Chip top layer thickness	2.0 mm						1.0 mm		1.0 mm		2.0 mm	
Chip base layer thickness	2.0 mm						3.0 mm		3.0 mm		2.0 mm	
Max. operating pressure	30 bar											
Material	B270 glass											
Fabrication process	HF etching and thermal bonding											
Channel coating	None (Hydrophilic)	Hydrophobic/Fluorophilic*	None (Hydrophilic)	Hydrophobic/Fluorophilic*	None (Hydrophilic)	Hydrophobic	None (Hydrophilic)	Hydrophobic/Fluorophilic*	None (Hydrophilic)	Hydrophobic/Fluorophilic*	None (Hydrophilic)	Hydrophobic
Compatible interface	H Interface (Part No. 3000155)						Top Interface 4-Way (4mm) (Part No. 3000109)					
Compatible connector	Linear Connector 4-Way (Part No. 3000024)											

\* The standard hydrophobic coating is not optimized for use with fluorocarbon oils. We suggest using a fluorophilic coating on the channel surface.

## 4 Geometry

### 4.1 3D Flow Focusing Chip - 65 $\mu\text{m}$ (P.N: 3200958, 3200959, 3200960)






Figure 2 External dimensions



### 4.2 3D Flow Focusing Chip - 100 $\mu\text{m}$ (P.N: 3200433, 3200434, 3200515)




Figure 3 External dimensions




Figure 4 Channel cross section





Figure 5 Junction geometry



*Figure 6 Nominal pore dimensions*

4.3 3D Flow Focusing Chip - 170  $\mu\text{m}$  (P.N: 3200835, 3200836, 3200837)




Figure 7 External dimensions




Figure 9 Channel cross section

Figure 8 Junction geometry

4.4 3D Flow Focusing Chip - 100  $\mu\text{m}$  - 2 Reagent (P.N: 3200435, 3200436)






Figure 10 External dimensions





*Figure 12 Junction geometry*



*Figure 13 Nominal pore dimensions*

4.5 3D Flow Focusing Chip – 30  $\mu\text{m}$  (P.N.: 3200934, 3200935, 3200936)

Figure 14 External dimensions





Figure 15 Junction geometry




Figure 16 Channel cross-section at the junction

4.6 3D Flow Focusing Chip - 14  $\mu\text{m}$  (P.N: 3200437, 3200438)




Figure 17 External dimensions




Figure 17 On-chip filter




Figure 18 Nominal pore dimensions

Figure 19 Junction geometry




Figure 181 Channel cross-section elsewhere






Figure 20 Channel cross-section at junction

## 5 Surface Coatings

The glass channel surface is naturally hydrophilic. This will form organic droplets in an aqueous carrier phase. To form aqueous droplets in an organic phase, hydrophobic or fluorophilic coatings are required.



The hydrophobic coating is resistant to organic solvents. The coating is vulnerable to acidic or basic solutions, such as even 0.1M NaOH.

A coating step in the fabrication differentiates between hydrophobic and fluorophilic chips. Depending on coating, users are able to create water-in-oil emulsion, oil-in-water emulsion, gas-in-water bubbles, or gas-in-oil bubbles.



## 6 Droplet Formation

The size, consistency, and production rate of droplet formation is a function of several physical parameters, including:

- Channel size
- Viscosity and surface tension of the various fluids
- Presence of surfactants
- Miscibility of the fluids
- Use of hydrophobic, fluorophilic or hydrophilic coating on channel walls
- Total flow rate and relative flow rate of each fluid
- Flow stability


The Mitos P-Pump (Part No. 3200016) has been designed to provide stable, pulse-free flow for droplet formation. In combination with a 3D Flow Focusing Chip it is possible to generate droplets with extremely consistent diameters (monodispersed).

The 100  $\mu\text{m}$  3D Flow Focusing Chip and the 2 Reagent 100  $\mu\text{m}$  3D Flow Focusing Chip typically form droplets of around 25 – 70 $\mu\text{m}$  in diameter. The 170  $\mu\text{m}$  3D Flow Focusing Chip forms droplets of around 50 – 140 $\mu\text{m}$  in diameter. The 30  $\mu\text{m}$  3D Flow Focusing Chip forms droplets of around 10 – 25  $\mu\text{m}$  in diameter. The 14  $\mu\text{m}$  3D Flow Focusing Chip forms droplets of around 5-12  $\mu\text{m}$  in diameter.

Chip Name/Part No.	Example image	Example droplet diameter
100 $\mu\text{m}$ 3D Flow Focusing Chip 3200433, 3200434, 3200515		65 $\mu\text{m}$
170 $\mu\text{m}$ 3D Flow Focusing Chip 3200835, 3200836, 3200837	Similar droplet formation as the 100 $\mu\text{m}$ 3D Flow Focusing	
2 Reagent 100 $\mu\text{m}$ 3D Flow Focusing Chip 3200435, 3200436		30 $\mu\text{m}$




30  $\mu\text{m}$  3D Flow Focusing Chip  
3200934/3200935/3200936




25  $\mu\text{m}$

14  $\mu\text{m}$  3D Flow Focusing Chip  
3200436/3200437



$\sim 5\mu\text{m}$



*Figure 19 70um droplets collected from the output of the 100um 3D Flow Focusing Chip*

It is possible to generate droplets larger than the maximum stated here, exceeding the channel depth. However, as the 3D Flow Focusing chip range is targeted at droplet fluids with channel-fouling properties, it is recommended in such cases to operate below the maximum sizes stated above to prevent droplets from coming into contact with the channel walls. When using more inert fluids, these limits need not apply.

Droplet production rates up to 12,000 droplets per second have been produced in a single channel. Dolomite has worked on many droplet forming projects, so please contact us for advice.

## 7 Custom Options

Dolomite can also offer additional customization, for example producing the same geometry in quartz. Other chip configurations are available on request. If you would like to generate droplets of a different size, operate with challenging fluids or perform complex droplet functions, Dolomite can design the junction geometry required. The range of Dolomite services available covers all aspects of the development process from characterization of liquids for droplet generation to the design of commercial instruments in the field of droplet microfluidics. Please contact Dolomite to discuss your application.

## 8 IP License

Dolomite is a licensee of Japan Science and Technology Agency (“JST”) under JST’s microdroplet generation technology. This enables our customers to purchase and use our droplet chips for R&D purposes without any restriction from this comprehensive IP family. Contact us for more information about licensing this IP for your custom application or chip design.